

## INSTRUCTION FOR USE

### ReGlass - Lc 2.5g

#### Light Curing Radiopaque Glass Ionomer Cavity Liner

#### COMPOSITION

Resin Modified Bio active Glass filler, Methacrylate esters, Photoinitiators and Paste Forming Agent

ReGlass - Lc is a light curing resin modified Bioactive Glass filler flowable composite for cavity linings, pit and fissure sealing and smaller surface lesions under ATR.

#### Properties and Method of Application

##### As Cavity Liner

ReGlass – Lc is inserted with applying needle cannula (recommended single use) directly into caries free prepared cavity area in 1-2 mm thickness, Photo polymerization of the material is conducted with ultra violet light of 400-500 nanometer wave length for 30-40 seconds depending upon depth of material. After setting of material becomes very hard so you can carefully remove and smooth base with rotary instruments for further preparation. Etching and bonding is advised depending upon type of your cavity and approximate distance with pulp tissue  
Do not use ReGlass – Lc on margin of cavity and enamel surface.

#### INDICATION

- As a Cavity liner under all restorative materials for better thermal insulators.
- In class V cavities for abrasion defects.

#### CONTRAINDICATION

Drug idiosyncrasy

Use with care if person is allergic to methacrylates and BHT

Do not use directly on exposed Pulp area as a temporary or permanent closure.

#### SIDE EFFECTS

No any known side effects as far as all terms of application, transportation and storage are observed.

#### PRESENTATION

ReGlass – Lc 2.5g syringe

Ref No. 005-04-250

#### STORAGE

Store at 5 °C to 25 °C, in a well sealed container Attention! Close materials tightly after use

The material must not be used after the expiry date Shelf life – 2.5 years

#### MANUFACTURER

SafeEndo Dental India Pvt Ltd

A/6/2, Kamdhenu Industrial Estate,

Opp. Gorwa Water tank, Gorwa, Vadodara-390016

E-mail: info@safeendo.org

www.safeendo.org

**Not For Medicinal use, For use in Dentistry By Dental Professional only**

005-04-250

Date of last revision : 2021-04

**SafeEndo**  
enhance your endo

+91 79-22139977 / 22163335   www.safeendo.org  info@safeendo.org  
www.safeendoindia.com

**SafeEndo**  
enhance your endo

+91 79-22139977 / 22163335   www.safeendo.org  info@safeendo.org  
www.safeendoindia.com