

INSTRUCTION FOR USE

IodoCure Paste

COMPOSITION

Iodoform, Calcium containing filler, Hydrophobic paste-forming component

PROPERTIES

Paste has a strong bactericidal effect.

Iodoform provides prolonged antimicrobial effect in root canal microflora. When in contact with tissues iodoform undergoes gradual decomposition evolving free iodine.

Hydrophobic paste-forming component provides paste's fluidity and plasticity, ensures hermetic sealing of root canal, protects soft tissues when inserting paste out of the apex.

Calcium containing filler stimulates hard tissues formation at the apex level, provides alkaline (pH 12.8) antimicrobial medium and gradual paste hardening under the influence of in-the-canal moisture.

Paste hardens in canal under the influence of in-the-canal moisture within 1-3 months. During hardening and after it the material can be easily removed out of the canal.

Being out of apex in soft tissues paste reabsorbs.

When filling of milk teeth paste reabsorbs as quickly as tooth resorption is held.

Paste can be applied with gutta-percha or polypropylene points.

IodoCure does not prevent from formation and developing of dental bud

In order to avoid tooth crown coloring do not fill whole pulp cavity with the paste

Paste is intended for use in aggressive biological medium with the temperature range between 32°C and 42°C (in oral cavity at the temperature (37±1)°C and humidity 60-90 %).

INDICATIONS

Paste is used as temporary filling of tooth root canals during staged treatment and as permanent filling of tooth root canals at the final stage of the treatment.

Paste **IodoCure** is used for :

- Treatment of acute and chronic forms of periodontitis
- Treatment of pulpitis, granulomatous periodontitis and granulating periodontitis
- For Apexification during treatment of tooth root canals with open apex in periodontitis
- Treatment of cystic lesion and chronic granuloma
- In case of root perforation and root resorption
- As permanent filling material for infected root canals after its treatment
- As permanent filling material for milk teeth

CONTRAINDICATIONS

Drug idiosyncrasy

Use with great care with patients allergic to iodine

Do not misuse

Warning! In case of allergy to iodine extract the material and stop the use of it

SIDE EFFECTS

As far as all terms of storing, transporting and application are observed, there are no side effects

METHODS OF USE

Put a tip on the syringe before use.

Prepare root canal using liquid for widening, antiseptic treatment and dehydrating liquids.

For treatment insert a tip into the canal and without too much pressure on syringe plunger fill the upper part of the canal with material, isolate with cotton swab and fill the cavity with temporary filling material (SmartTemp, SafeEndo.), without compression. Paste can stay in the cavity from 7-8 days to 1 month under care of the specialist. During follow up, in case of improvement, paste can be removed out of the canal. Clean the canal and using new tip, root needle or root canal filling instrument fill the canal with fresh portion of **IodoCure** paste, isolate with cavity liner and make the restoration of tooth crown.

For filling insert paste **IodoCure** into prepared canal, closer to apex, with the use of syringe with tip. After applying of cavity liner make the restoration of tooth crown.

For treatment of milk teeth remove carious dentine and necrotizing tissues first, then fill the canal with hardening paste **IodoCure**, isolate with temporary filling. During follow up extract upper part of temporary filling and fill the cavity

PRESENTATION

Paste (syringe) 2g/2x2g Ref No 004-02-202

STORAGE

Store at the temperature from 5°C to 25°C

Do not use after expiry date

Shelf life – 3 years

MANUFACTURER

SafeEndo Dental India Pvt Ltd

A/6/2, Kamdhenu Industrial Estate, Opp. Gorwa Water tank, Gorwa, Vadodara-390016

E-mail: info@safeendo.org

www.safeendo.org

Not For Medicinal use, For use in Dentistry By Dental Professional only

Ref No 004-02-202

Date of last revision: 2021-08

SafeEndo
enhance your endo



+91 79-22139977 / 22163335



www.safeendo.org
www.safeendoindia.com



info@safeendo.org

SafeEndo
enhance your endo



+91 79-22139977 / 22163335



www.safeendo.org
www.safeendoindia.com



info@safeendo.org