
BT96 Biological Indicator

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Quality certification Certificado de calidad Bionova[®] BT96

Plasma or Vaporized Hydrogen Peroxide sterilization / Esterilización por Plasma o Vapor de Peróxido de Hidrógeno *Geobacillus stearothermophilus* ATCC 7953

LOT

WH202

Heat shock Population / Población _____ **CFU / UFC**

D - value / Valor D _____ **sec./seg.** **(2.0 mg/L V_{H₂O₂}, 50 °C)**

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Survival time = (log₁₀ labeled population - 2) x labeled D-value

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Kill time = (log₁₀ labeled population + 4) x labeled D-value

System’s sensitivity / Sensibilidad del Sistema:

System’s sensitivity is determined as the difference between the number of positive indicators after 7-day incubation and false negative indicators (negative by fluorescence readout and visually positive) in relation to the number of 7-day positive indicators. / La sensibilidad del sistema se determina como la diferencia entre el número de indicadores positivos después de 7 días de incubación y los indicadores falsos negativos (negativos por lectura de fluorescencia y positivos visualmente) en relación al número de indicadores positivos luego de 7 días de incubación.

**(*) Sensibilidad /

(

N

p

Positives at 7 days

)

−
(

N

f

False negatives

)

X 100
≥
97%
Sensitivity

{\displaystyle \left({\frac {(N^{p}\ Positives at 7\ days)}{(N^{f}\ False negatives)}}{}}\right)X100\geq 97\%}**

Parameters determined at time of manufacture according to ISO 11138-1: 2017 and IRAM 37102-1: 1999 standards. The values shown are reproducible only under the same conditions under which they were determined.

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ISO and USP Compliant ATCC is a registered trademark of American Type Culture Collection	
<i>Terragene[®]</i> and <i>Bionova[®]</i> are registered trademarks of Terragene S.A. <i>STERRAD[®]</i> is a registered trademark of Advanced Sterilization Products. <i>V-PRO[®]</i> is a registered trademark of Steris Corporation. <i>STERIZONE[®]</i> is a registered trademark of TSO, Inc.	
 Lic. Adrián J. Rovetto Director Técnico Technical Director	
Use exclusivo para profesionales e Instituciones Sanitarias.	
Explanation of Symbols	 Manufacture Date
 Product designed for use with Vaporized Hydrogen Peroxide sterilization cycles	 Expiration Date
 Batch number	 Manufacturer
Intended Use Table	
Model	Vaporized Hydrogen Peroxide Cycles
	STERRAD 100S 54 minutos
BT96	STERRAD Standard and Express cycles
	V-Pro Max and Sterizone VP4
Fluorescence Read Time	pH Color Change
30 minutes	48 hours

EN Biological Indicator For plasma or vaporized Hydrogen Peroxide sterilization

Indications for Use
United States
Terragene[®] Bionova[®] SCBI (BT96) is a self-contained biological indicator inoculated with viable 10⁹ *Geobacillus stearothermophilus* spores and is intended for monitoring the efficacy of vaporized hydrogen peroxide sterilization processes. BT96 has Super Rapid readout at 30 minutes at 60 °C. See Intended Use Table for reference.

Outside the United States Terragene [®] Bionova [®] BT96 Fluorescence Super Rapid Readout Biological	

EN Biological Indicator For plasma or vaporized Hydrogen Peroxide sterilization

Indications for Use
United States
Terragene[®] Bionova[®] SCBI (BT96) is a self-contained biological indicator inoculated with viable 10⁹ *Geobacillus stearothermophilus* spores and is intended for monitoring the efficacy of vaporized hydrogen peroxide sterilization processes. BT96 has Super Rapid readout at 30 minutes at 60 °C. See Intended Use Table for reference.

Outside the United States
Terragene[®] Bionova[®] BT96 Fluorescence Super Rapid Readout Biological

Indicators have been designed for quick and easy monitoring of Plasma or Vaporized Hydrogen Peroxide sterilization cycles.

Device Description
Terragene[®] Bionova[®] BT96 Fluorescence Super Rapid Readout Biological Indicators are single-use Self-Contained Biological Indicators (SCBIs) that consist of a polypropylene tube, a spore carrier and a glass ampoule with a culture medium, enclosed with a colored cap. Each tube contains a population of *Geobacillus stearothermophilus* ATCC 7953 spores inoculated on a spore carrier, a plastic cap with holes and a barrier permeable to Plasma or Vaporized Hydrogen Peroxide. Each BT96 has a Process Indicator on label that changes from purple to green when exposed to hydrogen peroxide.

Precautions
WARNING: Do not use Terragene[®] Bionova[®] BT96 SCBI to control EO, Dry Heat, Formaldehyde or other sterilization processes different from the specified. Do not reuse SCBI.
WARNING: Place one or more SCBI in sterilizing hard-to-reach areas to ensure all areas of the chamber are sterilized. Please evaluate all load configurations to ensure ALL hard-to-reach areas have been identified, and place a SCBI in each of those locations.
WARNING: Do not reuse the sterilizer until the SCBI test result is negative.

Instructions for Use
1. Identify the Terragene[®] Bionova[®] BT96 SCBI by writing the sterilizer number (in case of having more than one), load number and processing date on the label.
2. Pack the SCBI along with materials to be sterilized in an appropriate package according to recommended sterilization practices. Place the package in those areas which are considered most inaccessible to the sterilizing agent (e.g., the center of the load and areas near the door).
3. Sterilize as usual.

4. After the sterilization process has finished, open the sterilizer door, wait five minutes and remove the SCBI from the package.
CAUTION: Wear safety glasses and gloves when removing the Terragene[®] Bionova[®] BT96 SCBI from the sterilized package.
WARNING: Do not crush or handle the SCBI excessively, since this might cause the glass ampoule to burst. Let the SCBI cool down until it reaches room temperature.

5. Check the Process Indicator on SCBI label. A color change to green indicates that the SCBI has been exposed to Hydrogen Peroxide.
IMPORTANT: This color change does not evidence the process effectiveness to achieve sterility. If the color of the Process Indicator has not changed, check the sterilization process.
6. Press the lid to seal the tube. Crush the ampoule contained in the SCBI with an individual ampoule crusher or with the ampoule crusher placed within the incubator’s incubation area. Then shake the tube down vigorously, with movements similar to those performed to lower the temperature in a mercury thermometer, until the medium reaches the base of the tube and soaks the spore carrier entirely. Finally, place the SCBI in the incubator.
IMPORTANT: Use a non-sterilized SCBI as a positive control at least once per day, when a sterilization cycle is run. The positive control ensures that correct incubation conditions were met, the capability of the medium to promote rapid growth, viability of spores has not been altered due to improper storage temperature, humidity or proximity to chemicals and proper functioning of Terragene[®] Bionova[®] Auto-Reader Incubators. Both, the positive control indicator and the processed indicator, should belong to the same batch.
7. Incubate the processed indicator and the positive control indicator in the appropriate Terragene[®] Bionova[®] Auto-Reader incubators for a maximum of 30 minutes at (60 ± 2) °C for super rapid readout.
NOTE: Holding time between sterilization and incubation should not exceed a 7-day period. Fluorescence detection by the reader (excitation 340-380 nm / emission 455-465 nm) means failure in the sterilization process. If no fluorescence is detected after 30 minutes of incubation, the result is negative. The positive control must give positive fluorescence readout. It is good practice to incubate a positive control for visual color change.

Record positive results and discard the SCBIs immediately, as it is indicated below.

Super Rapid Readout: 30 minutes
The Super Rapid Readout must be carried out in the appropriate Terragene[®] Bionova[®] Auto-Reader Incubators. Fluorescence is emitted after the reader stimulates the spore carrier with UV light. Final readout of negative results is readily available after 30 minutes of incubation. Fluorescence readout is an indirect measure of the viability of *Geobacillus stearothermophilus* spores after the sterilization process (positive results). Furthermore, a failure in the sterilization process can also become evident by culture medium color-change.
Dose response: The high sensitivity of the fluorescence results at 30 minutes, conventional incubation for color change of Terragene[®] Bionova[®] SCBI BT96 is not an advantage.

Visual confirmation: 48 hours
Optionally, you can perform a visual color-change confirmation after a 48-hour incubation. If the sterilization process has not been successful, culture medium will turn to yellow during incubation at (60 ± 2) °C, indicating the presence of living spores. If sterilization was successful, culture medium will remain purple after incubation. The positive control must show a color change from purple to yellow for results to be valid.

Readout time: 7 days
A 7-day readout is optional and may be routinely performed. This is an initial validation of the 30-minute readout. Fluorescence results are compared to the 7-day visual readouts.
NOTE: If 7-day readout is performed, a humidified environment will be required to prevent medium from drying out.

Monitoring frequency
Follow the policies and procedures with the monitoring frequency specified by the professional associations and/or standards corresponding to your country. As recommended practice, and to provide optimal patient safety, Terragene[®] recommends that each sterilization load be monitored with the appropriate biological indicator.
Storage
Store in a dark place at temperature between 10-30 °C and 30-80 % relative humidity. Do not freeze. Do not store biological indicators near sterilizing agents or other chemical products.
Shelf life
Biological Indicators have an expiration date of 2 years from the date of manufacture when stored at recommended conditions. Do not use indicators after their expiration date. Chemical Process Indicator on SCBI label have an expiration date of 2 years when used as part of SCBI.
Endpoint Stability Reaction: chemical indicator endpoint shall remain unchanged for a period of not less than 6 months when stored at previously indicated conditions.

Disposal
Discard biological indicators after use according to your country’s healthcare and safety regulations. The positive biological indicators can be autoclaved in a gravity air displacement steam sterilizer at 121 °C for 30 minutes, 132 °C for 15 minutes or 134 °C for 10 minutes; or in a dynamic air removal steam sterilizer at 132 °C for 4 minutes or 135 °C for 3 minutes.

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