CD29

Rev.11 / 01.2021

Multivariable Process Indicators

Result Reference Guide

| | |
|--|--|
| | Unprocessed |
| | Failures One or more sterilization parameters were not achieved. The materials cannot be considered sterile and must be processed again. The main failure causes are: <ul style="list-style-type: none"> Inadequate cycle duration Improper cycle selection Overloading / incorrect loading of the sterilizer Poor steam quality Sterilizer malfunction Blocking indicator with large metal mass |
| | Processed and Satisfactory |

NOTE: Reference colors exhibited in printed prospects as well as those displayed in our website and soft copies of other documents, do not necessarily represent the real color as shown in the actual indicators.

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Terragene® and Chemdye® are registered trademarks of Terragene S.A

EN Multivariable Process Indicators

For Steam sterilization processes

Indications for Use

United States
Terragene® Chemdye® CD29 is a Chemical Process Indicator intended for monitoring the efficacy of Steam sterilization processes. The chemical indicator changes from yellow to dark brown/black to indicate that the conditions of the cycle have been met. See Intended Use Table for reference.

Outside the United States

Chemdye® CD29 Chemical Indicator has been designed for monitoring Steam sterilization processes.

Device Description

Chemdye® CD29 Multivariable Process Indicators (Type 4 according to ISO 11140-1:2014 standard) are single-use chemical indicators that consist of one-side laminated punched paper double strips printed with indicator ink. These indicators have been designed to monitor Steam sterilization processes within loads, ensuring an adequate exposure to the sterilizing agent during the sterilization process and allowing to distinguish between processed and unprocessed items. Innovative highly sensitive yellow ink was developed to turn to dark brown/black when the indicator is exposed to Steam considering the critical process variables which provides additional quality assurance. Endpoint condition: 121 °C, 15 minutes; 134 °C, 3.5 minutes. Indicators 100 % free of toxic heavy metals.

Precautions

Do not expose Chemdye® CD29 Multivariable Process Indicators to Ethylene Oxide, Dry Heat or any sterilization process other than Steam. Do not reuse.

Instructions for Use

- Place a Chemdye® CD29 Multivariable Process Indicator within each package, bag or tray that will be sterilized by Steam. In rigid containers, place the indicator in every corner, or at least two diagonally opposite corners, to ensure that the sterilization agent has reached the most inaccessible areas.
- Perform the sterilization cycle.
- After the sterilization process has finished, remove the indicators from the sterilizer and analyze the results. For all the chemical indicators exposed to Steam, the ink must turn to the dark brown/black reference color when the sterilization parameters are reached. For reference color, please refer to Result Reference Guide.

Monitoring frequency

Follow facility Policies and Procedures which should specify a chemical indicator monitoring frequency compliant with professional association recommended practices and/or national guidelines and standards. As a good practice and to provide optimal patient safety, Terragene® recommends that every sterilization load be monitored with an appropriate chemical indicator.

Storage

Indicators should be stored in a dry place, protected from light at temperatures between 10-30 °C, 30-80 % relative humidity. Do not wet. Do not store near sterilizing agents.

Shelf Life

Chemical Indicators have an expiration date of 5 years from the date of manufacture when stored at recommended storage conditions. Both the expiration date and batch number of the product are specified on the package label and in each unit of the product. Do not use indicators after their expiration date. Endpoint Stability Reaction: chemical indicator endpoint shall remain unchanged for a period of 6 months when stored at previously indicated conditions.

Disposal

Discard chemical indicators after use along with paper waste, according to your country's healthcare and safety regulations.

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