# IT26-C

## Moving Front Chemical Integrator

For Steam sterilization processes.

### Usage

Control of Steam sterilization processes between 118 °C and 138 °C.

### Applicable regulation

Designed under Quality Management Systems standards ISO 13485:2016/NS-EN ISO 13485:2016.  
ISO 11140-1:2014

### Authorization

ANMAT (Argentinean National Administration of Drugs, Food and Medical Devices) PM 1614-4.

### Classification

Class 1, according to risk (ANMAT).

### FDA 510(k)

K191021

### Characteristics

Type 5 Integrating Indicator (according to ISO 11140-1).

Size: 52 mm x 21 mm.

Chemical Integrator consists of a paper wick and a steam and temperature sensitive chemical pellet contained in a paper/film/foil laminate. The chemical pellet melts and migrates as a dark bar along the paper wick. The migration is visible through a zone marked ACCEPT or REJECT, thus indicating whether sterilization conditions were met. The extent of migration depends on steam quality, time, and temperature. The ACCEPT result is reached when a theoretical spore population reaches its kill time, indicating integration condition has been reached. This condition is calibrated with the kill time of  $10^6$  *Geobacillus stearothermophilus* spore population, calculated in BIER (Biological Indicator Evaluator Resistometer).

Conditions: saturated steam at 121 °C, 128 °C, 135 °C.

Optionally, the readout of the chemical indicator can be done using Trazanto Automatic Quality Control and Traceability System (see instructions for use in device's manual).

### Environmental conditions during manufacture

T= 15-30 °C, RH = 30-80 %.

### Storage conditions

T= 10-30 °C, RH = 30-80 %, keep out of direct light.

### Transport conditions

Storage conditions should be strictly followed.

Transport in closed and reinforced boxes in order to avoid damages. The transport of this product does not represent a risk for health.

### Shelf life

5 years.

### Packing

100 or 250 units per aluminum foil bag.

### **Labelling**

On product's packing: product code and description, process for intended use, instructions for use, storage conditions, product presentation, classification according to regulation, manufacturer information, batch number, manufacture and expiration date, SV at 121 °C, 128 °C and 135 °C, barcode and datamatrix code.

### **Possible target markets**

Healthcare, Food, Pharmaceutical and Medical Industries.

### **Other important information**

Read product's instructions for use thoroughly before use.

### **Precautions**

Do not store the product near sterilizing agents, strong alkaline or acidic products such as cleaning/disinfecting agents.

Do not expose this product to EO, Dry Heat or any sterilization processes other than Steam.

Do not reuse the sterilizer until the dark bar has reached the ACCEPT zone.

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