# BT50

## Spore Strip Biological Indicator

For Steam and Formaldehyde sterilization processes.

### Usage

Monitoring 121 °C - 135 °C Steam and Formaldehyde Sterilization Processes.

### Applicable regulation

Designed under Quality Management System standards ISO 13485:2016/NS-EN ISO 13485:2016.

ISO 11138-1:2017, ISO 11138-3:2017 and ISO 11138-5:2017; IRAM 37102-1:1999 and IRAM 37102-3:1999.

### Authorization

ANMAT (Argentinean National Administration of Drugs, Food and Medical Devices) PM 1614-1.

### Classification

Class 1, according to risk (ANMAT).

### Characteristics

25.0 mm x 70.0 mm non-absorbent material package. Brown code.

5.0 mm x 25.0 mm filter paper strip.

$10^4$ ,  $10^5$ ,  $10^6$  or  $10^7$  *Geobacillus stearotherophilus* ATCC 7953 spores per strip.

Upon completion of the sterilization cycle, the spore strip should be tested in a suitable culture medium. We suggest the use of Bionova® MC20-2 or MC1020-2 culture media.

Incubation conditions:

24 hours, between 55-62 °C for Steam sterilization processes.

48 hours, between 55-62 °C for Formaldehyde sterilization processes.

For conventional culture media like TSB, incubation for 7 days between 55-62 °C is recommended.

For population  $\geq 10^5$  spores:

D-Value<sub>STEAM</sub>: not lower than 1.5 minutes. Conditions: 121 °C.

Z-Value<sub>STEAM</sub>: not lower than 6 °C.

D-Value<sub>FORM</sub>: not lower than 6.0 min. Conditions: 60 °C, 1 mol/l formaldehyde

### Environmental conditions during manufacture

T= 15-30 °C, RH = 30-80 %. Sterility conditions are necessary only during the inoculation process performed in laminar flow.

### Storage conditions

T= 10-30 °C, RH = 30-80 %.

### Transport conditions

Storage conditions should be strictly followed.

Products should be transported in closed and reinforced boxes in order to avoid damages.

The transport of this product does not represent any risk for human health.

### Shelf life

2 years.

**Packing**

100 units per aluminium foil bag.

**Labelling**

On product's envelope: product code and description, process for intended use, batch number, expiration date, bacterial strain and datamatrix code.

On product's packing: code and description of the product, process for intended use, strain and population of bacteria, batch number, manufacture and expiration date, presentation, regulation, storage conditions, manufacturer information, barcode and datamatrix code.

**Possible target markets**

Healthcare and Industry.

**Other important information**

Incubation at 60 °C in Bionova® IC10/20 incubator is recommended.  
Read product's instructions for use thoroughly before use.

**Precautions**

Do not store the product near sterilizing agents.

Do not expose this product to radiation, H<sub>2</sub>O<sub>2</sub> or any sterilization process other than Steam and/or Formaldehyde.

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