





Improve your digital workflow

One of the most accurate and fastest scanners on the market.

Free your hands

Automatic and adaptive calibration, manual calibration is not necessary any more.

Make your experience comfortable

Smaller and lighter while achieving a 16mm x 14mm high-definition field of view. More friendly to people with sensitive teeth. (*The air duct inside the body has been upgraded.)



Accuracy

As one of the most accurate scanners on the market, HELIOS 500 has reached the state-of-the-art accuracy level- the average deviation of a in-vivo full arch scan is around **20.0µm.**



Speed

Scan a full arch in less than 30 seconds with HELIOS 500.


Weight

Only **185 grams**, one of the most lightweight scanners on the market, well-balanced, ergonomically designed, and easy to carry.

Al technology



Al technology filters the soft and moving tissues, offering users a smooth scanning experience and generating a natural-looking treatment outcomes.



Plug&Scan





Smart Reminder

Automatic Software Update





Cloud Storage

Online 3D Viewer





Timely After-sales Service



Voice Note

Scanner Status Icons



Recommended System Requirements

Processor	Intel®CoreTM i7	12th Generation.	base frequency 2	2.4 GHz (or better)

Memory 16GB (or more) DDR4, frequency 2666 MHz (or better)

Hard Drive 256GB (or more) SSD

NVIDIA®GeForce®GTX 3050
Graphics Card

(or better)

15.6" FHD (1920 x 1080)

(or more)

Others USB 3.0 port

Operation system Windows 10/11

Optional Touch screen

cations

H≘LIOS 500

Accuracy	Full-arch	20.0µm	
	Scan area	16 x 14mm (L)	
		12 x 12mm (S)	
Handpiece	Dimensions	217x54x24mm	
	Weight	185g	
Tip	Tip size(front nort)	22x18mm(L)	
	Tip size(front part)	18.6x16.7mm(S)	
Connectivity	Length	1.8m/Handpiece detachable	
	Туре	USB 3.0	
Oalibustian	Accuracy calibration	Calibration-free	
Calibration	Color calibration	Calibration-free	
Data	Output format	STL/PLY	




Orikam Healthcare India Pvt Ltd

Gurugram Office:

3rd Floor, Arogyam Building, Subhash Nagar, New Railway Road, Gurugram–122001 (Hr.) Contact: +91 124 4361189, +91 98212 30028



info@orikam.in, sales@orikam.in



Mumbai Office:

OHI Marketing, 4143, Eaze Zone Mall, Malad (W), Opp. Rustomjee Ozone, Mumbai- 400064 (Mh.)

Contact: +91 98209 57308, +91 63649 87444



Scan to follow us on Facebook orikamhealthcareindia

